

thermo scientific



Chemicals

AcroSeal extra dry solvents

Your air- and moisture-sensitive solution

Find it at eu.fishersci.com

 **fisher scientific**
part of Thermo Fisher Scientific

The AcroSeal system: performance, convenience, safety

How to use AcroSeal packaging

AcroSeal™ packaging is easy to use. To get the best results, we recommend the following steps:

1. Puncture the septum
2. Withdraw the required amount
3. Replace the nitrogen blanket
4. Store upright to prevent prolonged product contact with the seal

Advanced tips:

1. Prepare the nitrogen balloon in advance and use two balloons, one inside the other
2. Purge the doubled balloons with nitrogen three times before use
3. Bend the needle to make it easier to use and control
4. Flush each needle and syringe with nitrogen three times before use
5. Use the reaction flask to flush the syringe and needle prior to use

To get the best results from the septum, use 18- to 21-gauge needles and puncture in a clockwise pattern to avoid piercing the septum in the same place.

If you require a large amount of solvent, remove the whole cap under an inert atmosphere and decant the desired amount of product from the bottle.



AcroSeal extra dry solvents

Choose from a wide range of extra dry solvents for moisture-sensitive and biochemistry applications. Deuterated solvents are also available in AcroSeal packaging, offering a low-cost alternative to ampule packaging and reducing the hazards associated with handling ampules.

Available in three different grades:

- Standard grade: Suitable for most applications
- Molecular sieves grade: Stored over molecular sieves for prolonged shelf life
- Supreme grade: Filtered through 0.2 micron filters for dust-free applications

Below is the complete range of Thermo Scientific extra-dry solvents in AcroSeal packaging.

IBS No.	Description	CAS No.	Package Size
10689082	Acetone, 99.8%, extra dry	67-64-1	1 L
10580772	Acetone, 99.8%, extra dry	67-64-1	100 mL
10203042	Acetonitrile, 99.9%, extra dry	75-05-8	1 L
10257473	Acetonitrile, 99.9%, extra dry	75-05-8	2.5 L
10193051	Acetonitrile, 99.9%, extra dry	75-05-8	100 mL
15101237	Acetonitrile, 99.9%, extra dry - FOR UCB ONLY -	75-05-8	500 mL
10222052	Acetonitrile, 99.9%, extra dry over molecular sieve	75-05-8	1 L
10174813	Acetonitrile, 99.9%, extra dry over molecular sieve	75-05-8	2.5 L
10353732	Acetonitrile, 99.9%, extra dry over molecular sieve	75-05-8	100 mL
10363413	Acetonitrile, 99.9%, extra dry over molecular sieve	75-05-8	500 mL
12317103	Acetonitrile, 99.9%, extra dry, package of 4 x 25 mL bottles	75-05-8	100 mL
11924721	Anisole, 99%, extra dry	100-66-3	1 L
11934721	Anisole, 99%, extra dry	100-66-3	100 mL
10721244	Benzyl alcohol, 98+%, extra dry	100-51-6	1 L
10731244	Benzyl alcohol, 98+%, extra dry	100-51-6	100 mL
10590374	Bis(2-methoxyethyl) ether, 99+%, extra dry	111-96-6	1 L
10500384	Bis(2-methoxyethyl) ether, 99+%, extra dry	111-96-6	100 mL
11954701	Bis(2-methoxyethyl) ether, 99+%, extra dry over molecular sieve	111-96-6	1 L

IBS No.	Description	CAS No.	Package Size
11964701	Bis(2-methoxyethyl) ether, 99+%, extra dry over molecular sieve	111-96-6	100 mL
10247893	1-Butanol, 99+%, extra dry	71-36-3	1 L
10483575	1-Butanol, 99+%, extra dry	71-36-3	100 mL
10653605	sec-Butanol, 99+%, extra dry	78-92-2	1 L
10537264	sec-Butanol, 99+%, extra dry	78-92-2	100 mL
10257563	2-Butanone, 99.5%, extra dry	78-93-3	1 L
10402895	2-Butanone, 99.5%, extra dry	78-93-3	100 mL
13107754	n-Butyl acetate, 99+%, extra dry	123-86-4	1 L
13177754	n-Butyl acetate, 99+%, extra dry	123-86-4	100 mL
10136043	tert-Butyl methyl ether, 99+%, extra dry	1634-04-4	1 L
13178453	tert-Butyl methyl ether, 99+%, extra dry	1634-04-4	100 mL
10089552	tert-Butyl methyl ether, 99+%, extra dry over molecular sieve	1634-04-4	1 L
10628273	tert-Butyl methyl ether, 99+%, extra dry over molecular sieve	1634-04-4	100 mL
10724124	Chlorobenzene, 99.8%, extra dry	108-90-7	1 L
10532574	Chlorobenzene, 99.8%, extra dry	108-90-7	100 mL
13197744	Chlorobenzene, 99.5%, extra dry over molecular sieve	108-90-7	1 L
13157784	Chlorobenzene, 99.5%, extra dry over molecular sieve	108-90-7	100 mL
13157794	1-Chlorobutane, 99.5%, anhydrous	109-69-3	1 L
13157714	1-Chlorobutane, 99.5%, anhydrous	109-69-3	100 mL
10139981	Chloroform, 99.9%, extra dry, stabilized	67-66-3	1 L
10399183	Chloroform, 99.9%, extra dry, stabilized	67-66-3	2.5 L
10071511	Chloroform, 99.9%, extra dry, stabilized	67-66-3	100 mL
10353582	Chloroform, 99.9%, extra dry over molecular sieve, stabilized	67-66-3	1 L
15322032	Chloroform, 99.9%, extra dry over molecular sieve, stabilized	67-66-3	2.5 L
10332702	Chloroform, 99.9%, extra dry over molecular sieve, stabilized	67-66-3	100 mL
15204228	Chloroform, 99.9%, extra dry over molecular sieve, stabilized	67-66-3	500 mL
10726305	Cyclohexane, 99.5%, extra dry	110-82-7	1 L
10550214	Cyclohexane, 99.5%, extra dry	110-82-7	2.5 L
10516061	Cyclohexane, 99.5%, extra dry	110-82-7	100 mL
10252572	Cyclohexane, 99.5%, extra dry over molecular sieve	110-82-7	1 L
10330474	Cyclohexane, 99.5%, extra dry over molecular sieve	110-82-7	2.5 L
10078891	Cyclohexane, 99.5%, extra dry over molecular sieve	110-82-7	100 mL
15244228	Cyclohexane, 99.5%, extra dry over molecular sieve	110-82-7	500 mL
13128493	Cyclopentane, 95+%, extra dry	287-92-3	1 L
13108563	Cyclopentane, 95+%, extra dry	287-92-3	100 mL
10453755	Cyclopentyl methyl ether, 99.5%, extra dry, stabilized	5614-37-9	1 L
10187243	Cyclopentyl methyl ether, 99.5%, extra dry, stabilized	5614-37-9	100 mL
11964461	Decahydronaphthalene, 99%, mixture of cis and trans, anhydrous	91-17-8	1 L
10791432	Decahydronaphthalene, 99%, mixture of cis and trans, anhydrous	91-17-8	100 mL
11934411	Decane, 99%, anhydrous	124-18-5	1 L
11944411	Decane, 99%, anhydrous	124-18-5	100 mL
10607764	Di-n-butyl ether, 99+%, extra dry	142-96-1	1 L
11934701	Di-n-butyl ether, 99+%, extra dry	142-96-1	2.5 L
10573324	Di-n-butyl ether, 99+%, extra dry	142-96-1	100 mL
10063443	1,2-Dichlorobenzene, 98+%, extra dry	95-50-1	1 L
10412895	1,2-Dichlorobenzene, 98+%, extra dry	95-50-1	100 mL
10612033	1,2-Dichloroethane, 99.8%, extra dry	107-06-2	1 L
11924471	1,2-Dichloroethane, 99.8%, extra dry	107-06-2	2.5 L
10333532	1,2-Dichloroethane, 99.8%, extra dry	107-06-2	100 mL
15342032	1,2-Dichloroethane, 99.5%, extra dry over molecular sieve	107-06-2	1 L
15352032	1,2-Dichloroethane, 99.5%, extra dry over molecular sieve	107-06-2	100 mL
16318227	Dichloromethane, anhydrous, 99.7+%, packaged under argon in resealable AcroSeal bottles, stabilized with amylene	75-09-2	250 mL
16328227	Dichloromethane, anhydrous, 99.7+%, packaged under argon in resealable AcroSeal bottles, stabilized with amylene	75-09-2	1 L
16338227	Dichloromethane, anhydrous, 99.7+%, packaged under argon in resealable AcroSeal bottles, stabilized with amylene	75-09-2	4 × 1 L
10477951	Dichloromethane, 99.9%, extra dry, stabilized	75-09-2	1 L
10437202	Dichloromethane, 99.9%, extra dry, stabilized	75-09-2	2.5 L
10429252	Dichloromethane, 99.9%, extra dry, stabilized	75-09-2	100 mL
10387841	Dichloromethane, 99.8%, extra dry over molecular sieve, stabilized	75-09-2	1 L
10295661	Dichloromethane, 99.8%, extra dry over molecular sieve, stabilized	75-09-2	2.5 L
10487532	Dichloromethane, 99.8%, extra dry over molecular sieve, stabilized	75-09-2	100 mL

IBS No.	Description	CAS No.	Package Size
10354263	Dichloromethane, 99.8%, extra dry over molecular sieve, stabilized	75-09-2	500 mL
12307083	Dichloromethane, 99.9%, extra dry, stabilized, package of 4 × 25 mL bottles	75-09-2	100 mL
11915161	Diethyl carbonate, 99%, anhydrous	105-58-8	1 L
11925161	Diethyl carbonate, 99%, anhydrous	105-58-8	100 mL
10169491	Diethyl ether, 99.5%, extra dry, stabilized	60-29-7	1 L
11974691	Diethyl ether, 99.5%, extra dry, stabilized	60-29-7	2.5 L
10018031	Diethyl ether, 99.5%, extra dry, stabilized	60-29-7	100 mL
10417372	Diethyl ether, 99.5%, extra dry over molecular sieve, stabilized	60-29-7	1 L
10501764	Diethyl ether, 99.5%, extra dry over molecular sieve, stabilized	60-29-7	2.5 L
10059031	Diethyl ether, 99.5%, extra dry over molecular sieve, stabilized	60-29-7	100 mL
10163893	Diethyl ether, 99.5%, extra dry over molecular sieve, stabilized	60-29-7	500 mL
12347103	Diethyl ether, 99.5%, extra dry, stabilized, package of 4 × 25 mL bottles	60-29-7	100 mL
10388333	N,N-Dimethylacetamide, 99.5%, extra dry over molecular sieve	127-19-5	1 L
10023743	N,N-Dimethylacetamide, 99.5%, extra dry over molecular sieve	127-19-5	100 mL
10424275	N,N-Dimethylacetamide, 99.5%, extra dry	127-19-5	1 L
13127804	N,N-Dimethylacetamide, 99.5%, extra dry	127-19-5	2.5 L
10699453	N,N-Dimethylacetamide, 99.5%, extra dry	127-19-5	100 mL
10483360	Dimethyl carbonate, 99+%, extra dry	616-38-6	1 L
13198473	Dimethyl carbonate, 99+%, extra dry	616-38-6	100 mL
16378227	N,N-Dimethylformamide, anhydrous, 99.8%, packaged under argon	68-12-2	4 × 100 ml
16388227	N,N-Dimethylformamide, anhydrous, 99.8%, packaged under argon	68-12-2	100 ml
16398227	N,N-Dimethylformamide, anhydrous, 99.8%, packaged under argon	68-12-2	1 L
16308237	N,N-Dimethylformamide, anhydrous, 99.8%, packaged under argon	68-12-2	4 × 1 L
10098721	N,N-Dimethylformamide, 99.8%, extra dry	68-12-2	1 L
10437392	N,N-Dimethylformamide, 99.8%, extra dry	68-12-2	2.5 L
10295761	N,N-Dimethylformamide, 99.8%, extra dry	68-12-2	100 mL
10534341	N,N-Dimethylformamide, 99.8%, extra dry over molecular sieve	68-12-2	1 L
10738742	N,N-Dimethylformamide, 99.8%, extra dry over molecular sieve	68-12-2	2.5 L
10045421	N,N-Dimethylformamide, 99.8%, extra dry over molecular sieve	68-12-2	100 mL
10088362	N,N-Dimethylformamide, 99.8%, extra dry over molecular sieve	68-12-2	500 mL
12337083	N,N-Dimethylformamide, 99.8%, extra dry, package of 4 × 25 mL bottles	68-12-2	100 mL
16368217	Dimethyl sulfoxide, HPLC grade, 99.9+%, packaged under argon	67-68-5	250 mL
16378217	Dimethyl sulfoxide, HPLC grade, 99.9+%, packaged under argon	67-68-5	1 L
16388217	Dimethyl sulfoxide, HPLC grade, 99.9+%, packaged under argon	67-68-5	4 × 1 L
16368207	1,4-Dioxane, anhydrous, 99.8%, packaged under argon	123-91-1	100 mL
16378207	1,4-Dioxane, anhydrous, 99.8%, packaged under argon	123-91-1	1 L
16388207	1,4-Dioxane, anhydrous, 99.8%, packaged under argon	123-91-1	4 × 1 L
10783203	1,4-Dioxane, 99.8%, extra dry, stabilized	123-91-1	1 L
10368143	1,4-Dioxane, 99.8%, extra dry, stabilized	123-91-1	2.5 L
10061931	1,4-Dioxane, 99.8%, extra dry, stabilized	123-91-1	100 mL
10352702	1,4-Dioxane, 99.5%, extra dry over Molecular Sieve	123-91-1	1 L
10388533	1,4-Dioxane, 99.5%, extra dry over Molecular Sieve	123-91-1	2.5 L
10762393	1,4-Dioxane, 99.5%, extra dry over Molecular Sieve	123-91-1	100 mL
15214228	1,4-Dioxane, 99.5%, extra dry over Molecular Sieve	123-91-1	500 mL
15249626	1,4-Dioxane, 99.8%, extra dry, stabilized, package of 4 × 25 mL bottles	123-91-1	100 mL
12140910	1,3-Dioxolane, 99.8%, anhydrous, stabilized with 75 ppm BHT	646-06-0	1 L
11986201	1,3-Dioxolane, 99.8%, anhydrous, stabilized with 75 ppm BHT	646-06-0	100 mL
11914411	n-Dodecane, 99%, anhydrous	112-40-3	1 L
11924411	n-Dodecane, 99%, anhydrous	112-40-3	100 mL
10483195	Ethanol, 99.5%, extra dry, absolute	64-17-5	1 L
10646134	Ethanol, 99.5%, extra dry, absolute	64-17-5	2.5 L
10617864	Ethanol, 99.5%, extra dry, absolute	64-17-5	100 mL
10736305	Ethyl acetate, 99.9%, extra dry	141-78-6	1 L
10698054	Ethyl acetate, 99.9%, extra dry	141-78-6	2.5 L
10469252	Ethyl acetate, 99.9%, extra dry	141-78-6	100 mL
11396515	Ethyl acetate, 99.9%, extra dry over molecular sieve	141-78-6	1 L
10732833	Ethyl acetate, 99.9%, extra dry over molecular sieve	141-78-6	2.5 L
10394732	Ethyl acetate, 99.9%, extra dry over molecular sieve	141-78-6	100 mL
11914451	Ethylbenzene, 99.8%, anhydrous	100-41-4	1 L
13177724	Ethylbenzene, 99.8%, anhydrous	100-41-4	100 mL

IBS No.	Description	CAS No.	Package Size
11914631	Ethylene glycol, 99.8%, anhydrous	107-21-1	1 L
11934631	Ethylene glycol, 99.8%, anhydrous	107-21-1	100 mL
11994721	Ethylene glycol dimethyl ether, 99.5%, extra dry over molecular sieve	110-71-4	1 L
11904731	Ethylene glycol dimethyl ether, 99.5%, extra dry over molecular sieve	110-71-4	100 mL
10614375	Ethylene glycol dimethyl ether, 99+%, extra dry	110-71-4	1 L
10336744	Ethylene glycol dimethyl ether, 99+%, extra dry	110-71-4	100 mL
11954891	Ethyl formate, 98+%, extra dry	109-94-4	100 mL
10580994	n-Heptane, 99+%, extra dry	142-82-5	1 L
10396591	n-Heptane, 99+%, extra dry	142-82-5	100 mL
10500542	n-Heptane, 99+%, extra dry over molecular sieve	142-82-5	1 L
11326525	n-Heptane, 99+%, extra dry over molecular sieve	142-82-5	2.5 L
10783403	n-Heptane, 99+%, extra dry over molecular sieve	142-82-5	100 mL
11954411	n-Hexadecane, 99%, anhydrous	544-76-3	1 L
11964411	n-Hexadecane, 99%, anhydrous	544-76-3	100 mL
10292412	n-Hexane, 96+%, extra dry	110-54-3	1 L
10359233	n-Hexane, 96+%, extra dry	110-54-3	2.5 L
10061871	n-Hexane, 96+%, extra dry	110-54-3	100 mL
10689272	n-Hexane, 97%, extra dry over molecular sieve	110-54-3	1 L
15332032	n-Hexane, 97%, extra dry over molecular sieve	110-54-3	2.5 L
10437402	n-Hexane, 97%, extra dry over molecular sieve	110-54-3	100 mL
15224228	n-Hexane, 97%, extra dry over molecular sieve	110-54-3	500 mL
10431484	Hexanes, 99.9%, extra dry, mixture of isomers	92112-69-1	1 L
10401674	Hexanes, 99.9%, extra dry, mixture of isomers	92112-69-1	100 mL
11974611	Hexyl alcohol, 99%, anhydrous	111-27-3	1 L
11984611	Hexyl alcohol, 99%, anhydrous	111-27-3	100 mL
10058701	Isopropanol, 99.8%, extra dry	67-63-0	1 L
10533704	Isopropanol, 99.8%, extra dry	67-63-0	2.5 L
10787962	Isopropanol, 99.8%, extra dry	67-63-0	100 mL
10743793	Isopropanol, 99.5%, extra dry over molecular sieve	67-63-0	1 L
10227163	Isopropanol, 99.5%, extra dry over molecular sieve	67-63-0	2.5 L
10364162	Isopropanol, 99.5%, extra dry over molecular sieve	67-63-0	100 mL
15881022	Isopropanol, 99.5%, extra dry over molecular sieve	67-63-0	500 mL
16398217	Methanol, anhydrous, 99.9%, packaged under argon	67-56-1	250 mL
16308227	Methanol, anhydrous, 99.9%, packaged under argon	67-56-1	1 L
10511732	Methanol, 99.9%, extra dry	67-56-1	1 L
11984591	Methanol, 99.9%, extra dry	67-56-1	2.5 L
10747582	Methanol, 99.9%, extra dry	67-56-1	100 mL
10649492	Methanol, 99.8%, extra dry over molecular sieve	67-56-1	1 L
10468283	Methanol, 99.8%, extra dry over molecular sieve	67-56-1	2.5 L
10323442	Methanol, 99.8%, extra dry over molecular sieve	67-56-1	100 mL
11366525	Methanol, 99.8%, extra dry over molecular sieve	67-56-1	500 mL
12337103	Methanol, 99.9%, extra dry, package of 4 x 25 mL bottles	67-56-1	100 mL
10793154	2-Methoxyethanol, 99+%, extra dry	109-86-4	1 L
13118633	2-Methoxyethanol, 99+%, extra dry	109-86-4	100 mL
10195512	Methyl acetate, 99+%, extra dry	79-20-9	1 L
10759873	Methyl acetate, 99+%, extra dry	79-20-9	2.5 L
10037772	Methyl acetate, 99+%, extra dry	79-20-9	100 mL
15711819	2-Methylbutane, 99+%, extra dry	78-78-4	1 L
13111994	2-Methylbutane, 99+%, extra dry	78-78-4	100 mL
13171914	Methylcyclohexane, 98+%, extra dry	108-87-2	1 L
13172024	Methylcyclohexane, 98+%, extra dry	108-87-2	100 mL
13294663	2-Methyl-1-propanol, 99+%, extra dry	78-83-1	1 L
13128693	2-Methyl-1-propanol, 99+%, extra dry	78-83-1	100 mL
10466622	1-Methyl-2-pyrrolidinone, 99.5%, extra dry	872-50-4	1 L
10533864	1-Methyl-2-pyrrolidinone, 99.5%, extra dry	872-50-4	2.5 L
10229441	1-Methyl-2-pyrrolidinone, 99.5%, extra dry	872-50-4	100 mL
10079051	1-Methyl-2-pyrrolidinone, 99.5%, extra dry over molecular sieve	872-50-4	1 L
10304122	1-Methyl-2-pyrrolidinone, 99.5%, extra dry over molecular sieve	872-50-4	2.5 L
10383692	1-Methyl-2-pyrrolidinone, 99.5%, extra dry over molecular sieve	872-50-4	100 mL
15234228	1-Methyl-2-pyrrolidinone, 99.5%, extra dry over molecular sieve	872-50-4	500 mL
10774503	Methyl sulfoxide, 99.7+%, extra dry	67-68-5	1 L
10308103	Methyl sulfoxide, 99.7+%, extra dry	67-68-5	2.5 L

IBS No.	Description	CAS No.	Package Size
10081951	Methyl sulfoxide, 99.7+%, extra dry	67-68-5	100 mL
10225631	Methyl sulfoxide, 99.7+%, extra dry over molecular sieve	67-68-5	1 L
10545111	Methyl sulfoxide, 99.7+%, extra dry over molecular sieve	67-68-5	2.5 L
10487172	Methyl sulfoxide, 99.7+%, extra dry over molecular sieve	67-68-5	100 mL
10659643	Methyl sulfoxide, 99.7+%, extra dry over molecular sieve	67-68-5	500 mL
10409193	2-Methyltetrahydrofuran, 99+%, extra dry, stabilized	96-47-9	1 L
11976171	2-Methyltetrahydrofuran, 99+%, extra dry, stabilized	96-47-9	2.5 L
10114003	2-Methyltetrahydrofuran, 99+%, extra dry, stabilized	96-47-9	100 mL
10741354	2-Methyltetrahydrofuran, 99+%, extra dry over molecular sieve, stabilized	96-47-9	1 L
10104243	2-Methyltetrahydrofuran, 99+%, extra dry over molecular sieve, stabilized	96-47-9	100 mL
10676314	2-Methyltetrahydrofuran, 99+%, extra dry, stabilizer free	96-47-9	1 L
10550394	2-Methyltetrahydrofuran, 99+%, extra dry, stabilizer free	96-47-9	100 mL
10755642	n-Nonane, 99%, anhydrous	111-84-2	100 mL
10592374	n-Octane, 99+%, extra dry	111-65-9	1 L
13188633	n-Octane, 99+%, extra dry	111-65-9	100 mL
11984601	1-Octanol, 99%, anhydrous	111-87-5	1 L
13162164	1-Octanol, 99%, anhydrous	111-87-5	100 mL
10443185	n-Pentane, 99+%, extra dry	109-66-0	1 L
10279333	n-Pentane, 99+%, extra dry	109-66-0	100 mL
11974411	n-Pentane, 99+%, extra dry over molecular sieve	109-66-0	1 L
11984411	n-Pentane, 99+%, extra dry over molecular sieve	109-66-0	100 mL
10607954	Petroleum ether, extra dry, boiling range 40 to 60°C, water <50ppm	64742-49-0	1 L
13178643	Petroleum ether, extra dry, boiling range 40 to 60°C, water <50ppm	64742-49-0	100 mL
10592944	1-Propanol, 99.5%, extra dry	71-23-8	1 L
10003493	1-Propanol, 99.5%, extra dry	71-23-8	100 mL
13118663	Propylene carbonate, 99.5%, anhydrous	108-32-7	100 mL
13177714	Propylene carbonate, 99.5%, anhydrous	108-32-7	800 mL
16398207	Pyridine, anhydrous, 99.5+%, packaged under argon	110-86-1	100 mL
16308217	Pyridine, anhydrous, 99.5+%, packaged under argon	110-86-1	1 L
10502192	Pyridine, 99.5%, extra dry	110-86-1	1 L
10400374	Pyridine, 99.5%, extra dry	110-86-1	2.5 L
10113641	Pyridine, 99.5%, extra dry	110-86-1	100 mL
10353092	Pyridine, 99.5%, extra dry over molecular sieve	110-86-1	1 L
10781974	Pyridine, 99.5%, extra dry over molecular sieve	110-86-1	2.5 L
10700693	Pyridine, 99.5%, extra dry over molecular sieve	110-86-1	100 mL
12834741	Pyridine, 99.5%, extra dry over molecular sieve	110-86-1	500 mL
10168751	Tetrahydrofuran, 99.85%, extra dry, stabilized	109-99-9	1 L
10502232	Tetrahydrofuran, 99.85%, extra dry, stabilized	109-99-9	2.5 L
10613372	Tetrahydrofuran, 99.85%, extra dry, stabilized	109-99-9	100 mL
10798552	Tetrahydrofuran, 99.5%, extra dry over molecular sieve, stabilized	109-99-9	1 L
10256061	Tetrahydrofuran, 99.5%, extra dry over molecular sieve, stabilized	109-99-9	2.5 L
10292182	Tetrahydrofuran, 99.5%, extra dry over molecular sieve, stabilized	109-99-9	100 mL
10641304	Tetrahydrofuran, 99.5%, extra dry over molecular sieve, stabilized	109-99-9	500 mL
12317083	Tetrahydrofuran, 99.85%, extra dry, stabilized, package of 4 × 25 mL bottles	109-99-9	100 mL
15245508	Tetrahydrofuran, 99.9%, extra dry, unstabilized	109-99-9	1 L
15370627	Tetrahydrofuran, 99.9%, extra dry, unstabilized	109-99-9	2.5 L
15265508	Tetrahydrofuran, 99.9%, extra dry, unstabilized	109-99-9	100 mL
16318217	Tetrahydrofuran, anhydrous, 99.8+%, BHT-free, over molecular sieves, packed with argon	109-99-9	100 mL
16328217	Tetrahydrofuran, anhydrous, 99.8+%, BHT-free, over molecular sieves, packed with argon	109-99-9	1 L
16348227	Tetrahydrofuran, anhydrous, 99.8+%, Stabilized with 0.025% BHT, packaged under argon	109-99-9	100 mL
16358227	Tetrahydrofuran, anhydrous, 99.8+%, Stabilized with 0.025% BHT, packaged under argon	109-99-9	1 L
16368227	Tetrahydrofuran, anhydrous, 99.8+%, Stabilized with 0.025% BHT, packaged under argon	109-99-9	4 × 1 L
16338217	Tetrahydrofuran, spectrophotometric grade, 99.7+%, unstabilized, packaged under argon	109-99-9	250 ml
16348217	Tetrahydrofuran, spectrophotometric grade, 99.7+%, unstabilized, packaged under argon	109-99-9	1 L

IBS No.	Description	CAS No.	Package Size
16358217	Tetrahydrofuran, spectrophotometric grade, 99.7+%, unstabilized, packaged under argon	109-99-9	4 × 1 L
11974461	1,2,3,4-Tetrahydronaphthalene, 99%, extra dry	119-64-2	1 L
11984461	1,2,3,4-Tetrahydronaphthalene, 99%, extra dry	119-64-2	100 mL
10109731	Toluene, 99.85%, extra dry	108-88-3	1 L
10419633	Toluene, 99.85%, extra dry	108-88-3	2.5 L
10717202	Toluene, 99.85%, extra dry	108-88-3	100 mL
10510872	Toluene, 99.85%, extra dry over molecular sieve	108-88-3	1 L
10150662	Toluene, 99.85%, extra dry over molecular sieve	108-88-3	2.5 L
10477172	Toluene, 99.85%, extra dry over molecular sieve	108-88-3	100 mL
10354113	Toluene, 99.85%, extra dry over molecular sieve	108-88-3	500 mL
12327103	Toluene, 99.85%, extra dry, package of 4 × 25 mL bottles	108-88-3	100 mL
11914481	Trichloroethylene, 99%, extra dry, stabilized	79-01-6	1 L
11924481	Trichloroethylene, 99%, extra dry, stabilized	79-01-6	100 mL
13102174	Triethyl orthoformate, 98%, anhydrous	122-51-0	100 mL
11974941	Triethyl orthoformate, 98%, anhydrous	122-51-0	800 mL
11984941	Trimethyl orthoformate, 99+%, anhydrous	149-73-5	100 mL
11994941	Trimethyl orthoformate, 99+%, anhydrous	149-73-5	800 mL
10565714	2,2,4-Trimethylpentane, 99.5%, extra dry	540-84-1	1 L
10732823	2,2,4-Trimethylpentane, 99.5%, extra dry	540-84-1	100 mL
10692613	2,2,4-Trimethylpentane, 99.5%, extra dry over molecular sieve	540-84-1	1 L
10078901	2,2,4-Trimethylpentane, 99.5%, extra dry over molecular sieve	540-84-1	100 mL
11904441	o-Xylene, 99%, extra dry	95-47-6	1 L
11914441	o-Xylene, 99%, extra dry	95-47-6	100 mL
11924441	m-Xylene, 99%, extra dry	108-38-3	1 L
11944441	m-Xylene, 99%, extra dry	108-38-3	100 mL
10453465	Xylenes, 98+%, extra dry, mixed isomers	1330-20-7	1 L
13188473	Xylenes, 98+%, extra dry, mixed isomers	1330-20-7	100 mL

Browse the complete portfolio of products supporting research and manufacturing.

For more information visit: eu.fishersci.com/go/thermo-scientific-chemicals

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at **Belgium:** fishersci.be **Denmark:** fishersci.dk
Germany: fishersci.de **Ireland:** fishersci.ie **Italy:** fishersci.it
Finland: fishersci.fi **France:** fishersci.fr **Netherlands:** fishersci.nl
Norway: fishersci.no **Portugal:** fishersci.pt **Spain:** fishersci.es
Sweden: fishersci.se **Switzerland:** fishersci.ch **UK:** fishersci.co.uk

© 2021 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.