



Fisher Chemical Aqualine Reagents

Your solution for water content
determination by Karl Fischer titration

COULOMETRIC TITRATION
VOLUMETRIC TITRATION
WATER STANDARD

AQUALINE REAGENTS

Our Fisher Chemical™ Aqualine™ Karl Fischer reagents portfolio is designed to meet the needs of the analytical chemist by providing accurate water content determination using volumetric or coulometric titration with unique benefits.

AQUALINE COULOMETRIC RANGE: FOR LOW WATER CONTENT AT PPM LEVEL

Fisher Chemical Aqualine coulometric reagents are ideal for use in coulometric Karl Fischer titrations for detecting low concentrations of water. Our Aqualine anolyte and catholyte solutions have been re-formulated to offer better performance. The improved formulation increases both the speed and accuracy of titration when determining water content at microgram level.

Highlights:

- **Fast** – Reach the endpoint quickly
- **Convenient** – Long product shelf-life
- **Reliable** – Very stable endpoint



AQUALINE COULOMETRIC RANGE PERFORMANCE

The performance of Fisher Chemical Aqualine Electrolyte A coulometric reagent was tested in terms of water recovery and titration. Results indicate that Aqualine coulometric reagents are fast and accurate with a stable endpoint.




Figure 1: The speed of titration of a 1mL injection of methanol with Aqualine Electrolyte A was measured.

AQUALINE WATER STANDARD

We offer a series of long shelf-life Aqualine water standard reagents, which support the calibration of the Karl Fischer titrator instrument. Our water standard reagents are packaged in glass ampoules for your convenience.





AQUALINE VOLUMETRIC RANGE: FOR HIGH WATER CONTENT ANALYSIS

For the Karl Fischer titration by volumetry, we offer you the choice between single component reagents, two component reagents and reagents especially designed for aldehydes and ketones. **Our Aqualine reagents for Karl Fischer titration by volumetry have comparable performance to other Karl Fischer reagents on the market.**

Highlights:

- **Reliable** – Fast and stable endpoints ensure reliable and accurate results
- **Safe** – Low toxicity and pyridine free
- **Convenient** – Available as a one or two component solution

AQUALINE VOLUMETRIC RANGE PERFORMANCE

Fisher Chemical Aqualine Complete 5 water recovery and speed titration performance were tested, using imidazole as base.

Speed titration performance




Figure 3: The speed of titration of 50µg water with Aqualine Complete 5 was recorded using an automated volumetric moisture meter.

RESULTS OF TESTS ON ACTUAL SAMPLES

Sample	Working medium	Aqualine	
		Titration time min (std dev)	Water content % (std dev)
Acetic acid (water added)	Methanol	3:23 (0:12)	1.51 (0.00)
Triethylamine	Methanol	1:24 (0:17)	1.39 (0.01)
Dichloroacetic acid	Methanol + Imidazole	1:11 (0:13)	1.90 (0.01)
Vegetable oil	50% Methanol + 50% Xylene	1:39 (0:28)	0.03 (0.00)
Shower gel	Methanol	2:01 (0:02)	81.09 (0.06)
Acetone	Ketosolver	1:50 (0:13)	1.01 (0.01)
Hand cream	50% Methanol + 50% Chloroform	1:24 (0:08)	84.80 (0.03)
Coffee	50% Methanol + 50% Formamide	5:30 (0:01)	5.44 (0.01)
Chocolate	50% Methanol + 50% Chloroform	4:26 (0:31)	1.75 (0.02)

SELECT THE SUITABLE FISHER CHEMICAL AQUALINE REAGENT FOR YOUR KARL FISCHER TITRATION

Aqualine coulometric range: for low water content at ppm level

Product code	Product Description	Manufacturer product code	Pack size
Anolyte solutions			
These contain methanol and chloroform as solvents			
10302912	Aqualine Electrolyte A – For general use in conventional cells with a diaphragm	K/2500/08	500mL
10222132	Aqualine Electrolyte AD – For use in fritless (diaphragm-free) cells	K/2510/08	500mL
These contain methanol as a solvent			
10689672	Aqualine Electrolyte AD-G – For use in fritless (diaphragm-free) cells	K/2515/08	500mL
10724893	Aqualine Electrolyte AG – For general use in conventional cells with diaphragm	K/2520/08	500mL
This contains methanol and pentan-1-ol as solvents			
10592785	Aqualine Electrolyte AG-H – For samples with high level of hydrocarbon content	K/2530/08	500mL
Catholyte solution			
10169371	Aqualine Electrolyte CG – For general use in conventional cells with a diaphragm, contains methanol as a solvent	K/2560/04	25mL

Aqualine water standard

10609103	Aqualine Standard 0.2 – 0.2 mg/ml H ₂ O standard	K/2740/99	10 x 4mL
10323102	Aqualine Standard 1.0 – 1 mg/ml H ₂ O standard	K/2710/99	10 x 4mL
10609672	Aqualine Standard 5.0 – 5 mg/ml H ₂ O standard	K/2730/08	500mL
10546242	Aqualine Standard 10.0 – 10 mg/ml H ₂ O standard	K/2720/99	10 x 8mL
10723383	Aqualine Water Standard-KF oven – For use in KF oven technique, contains 5.55 ±0.05% water	K/2760/45	10g
10261552	Aqualine Sodium-tartrate Dihydrate Standard – Primary standard for volumetric analysis, contains 15.66 ±0.05% water	K/2770/48	100g

Aqualine volumetric range: for high water content analysis

Single component reagents			
10538810	Aqualine Complete 1 – Water equivalent 1mg H ₂ O/ml	K/1900/15	1L
10744761	Aqualine Complete 2 – Water equivalent 2mg H ₂ O/ml	K/1950/15	1L
10754761	Aqualine Complete 2 – Water equivalent 2mg H ₂ O/ml	K/1950/17	2.5L
10181570	Aqualine Complete 5 – Water equivalent 5mg H ₂ O/ml	K/2000/15	1L
10676262	Aqualine Complete 5 – Water equivalent 5mg H ₂ O/ml	K/2000/17	2.5L
Reagents for aldehydes and ketones			
10092042	Aqualine Complete 5K – Water equivalent 5 mg H ₂ O/ml	K/2250R/15	1L
10205592	Aqualine Matrix K – Matrix K should be used in conjunction with Complete 5K	K/2300R/15	1L
Two component reagents			
10786291	Aqualine Titrant 2 – Water equivalent 2 mg H ₂ O/ml	K/2150/15	1L
10284100	Aqualine Titrant 2 – Water equivalent 2 mg H ₂ O/ml	K/2150/17	2.5L
10528810	Aqualine Titrant 5 – Water equivalent 5 mg H ₂ O/ml	K/2200/15	1L
10172050	Aqualine Titrant 5 – Water equivalent 5 mg H ₂ O/ml	K/2200/17	2.5L
10101580	Aqualine Solvent	K/2100/15	1L
10264390	Aqualine Solvent	K/2100/17	2.5L
10199511	Aqualine Solvent CM – Solvent for samples with high hydrocarbon content	K/2110/15	1L
10510852	Aqualine Solvent CM – Solvent for samples with high hydrocarbon content	K/2110/17	2.5L

All our Aqualine solutions are pyridine free

To place an order, contact your local Fisher Scientific office.



©2016 Thermo Fisher Scientific Inc. All rights reserved.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

EUROPE

Austria
info.austria@thermofisher.com
www.at.fishersci.com
Tel: 0800 20 88 40
Fax: 0800 20 66 90

Belgium
be.fisher@thermofisher.com
www.be.fishersci.com
Tel: 056 260 260
Fax: 056 260 270

Czech Republic
info.cz@thermofisher.com
www.cz.fishersci.com
Tel: 466 798 230
Fax: 466 435 008

Denmark
kundeservice@thermofisher.com
www.dk.fishersci.com
Tel: 70 27 99 20
Fax: 70 27 99 29

Finland
fisher.fi@thermofisher.com
www.fi.fishersci.com
Tel: 09-802 76 280
Fax: 09-802 76 235

France
fr.fisher@thermofisher.com
www.fr.fishersci.com
Tel: 03 88 67 14 14
Fax: 03 88 67 11 68

Germany
info.germany@thermofisher.com
www.de.fishersci.com
Tel: 02304 932-5
Fax: 02304 932-950

Ireland
fsie.sales@thermofisher.com
www.ie.fishersci.com
Tel: 01 885 5854
Fax: 01 899 1855

Italy
it.fisher@thermofisher.com
www.it.fishersci.com
Tel: 02 950 59 478
Fax: 02 950 59 479

The Netherlands
nl.info@thermofisher.com
www.nl.fishersci.com
Tel: 020 487 70 00
Fax: 020 487 70 70

Norway
fisher.no@thermofisher.com
www.no.fishersci.com
Tel: 22 95 59 59
Fax: 22 95 59 40

Portugal
pt.fisher@thermofisher.com
www.pt.fishersci.com
Tel: 21 425 33 50
Fax: 21 425 33 51

Spain
es.fisher@thermofisher.com
www.es.fishersci.com
Tel: 902 239 303
Fax: 902 239 404

Sweden
fisher.se@thermofisher.com
www.se.fishersci.com
Tel: 31 352 32 00
Fax: 31 352 32 50

Switzerland
info.ch@thermofisher.com
www.ch.fishersci.com
Tel: 056 618 41 11
Fax: 056 618 41 41

United Kingdom
fsuk.sales@thermofisher.com
www.fisher.co.uk
Tel: 01509 555500
Fax: 01509 555111

ASIA

India
qfc.customercare@thermofisher.com
www.fishersci.in
Tel: +91-22-6680 3000
Fax: +91-22-6680 3001

Korea
sales.fsk@thermofisher.com
www.kr.fishersci.com
Tel: 02-3420-8700
Fax: 02-3420-8710

Malaysia
enquiry.my@thermofisher.com
www.fishersci.com.my
Tel: (603) 51228888
Fax: (603) 5121 8899

Singapore
enquiry.sg@thermofisher.com
www.fishersci.com.sg
Tel: (65) 6873 6006
Fax: (65) 6873 5005

OCEANIA

Australia
AUinfo@thermofisher.com
www.thermofisher.com.au
Tel: 1300-735-292
Fax: 1800-067-639

New Zealand
NZinfo@thermofisher.com
www.thermofisher.co.nz
Tel: 0800-933-966
Fax: +64 9 980 6788



A Thermo Fisher Scientific Brand